

# AnteAGE® Home Microneedling Kit

## Roll Your Way to Healthy Skin



Achieve incredible results with this natural approach to radiant skin. Stem cell growth factors maximize the regenerative power of microneedling while high grade hyaluronic acid locks in hydration. AnteAGE® brings it all together



### Instructions:

Disinfect roller by soaking in rubbing alcohol before and after each use.



Apply AnteAGE® Home Microneedling Solution to treatment area.



Take the roller and roll back and forth over the treatment area. Do this horizontally, vertically, and diagonally.



Apply AnteAGE® Home Microneedling Solution to treatment area again once rolling is complete.

### Solution Ingredients:

Water (Aqua), Human Bone Marrow Stem Cell Conditioned Media, Hyaluronic Acid, Benzyl Alcohol, Dehydroacetic Acid